

# Giulia

## Giulia

### Neonatal Ventilator



# Giulia

## A new challenge in neonatal respiratory care

In the multifactorial pathogenesis of Broncho-Pulmonary Dysplasia (BPD), ventilator induced lung-injury (VILI) is considered to be a significant factor. This prompted the development of new non-invasive respiratory techniques, which would be more effective than nasal CPAP, for the treatment of the newborn with RDS. Nasal intermittent positive pressure ventilation is a non-invasive mode of ventilation that combines nasal CPAP with some intermittent mandatory breaths. It may be non-synchronized (NIPPV), or synchronized (SNIPPV) with the infant's breathing efforts. At GINEVRI we have developed a revolutionary flow-sensor to carry out SNIPPV. Figure 1 shows this flow sensor being used for the treatment of a baby of 650 grams: the device is very reliable, comfortable and easy to fit.

Several clinical trials favor SNIPPV, probably because delivering the inspiratory pressure immediately after the start of a respiratory effort, when the glottis is open, allows pressure to be transmitted effectively to the lungs. Figures 2 and 3 show the different ways that NIPPV and flow-SNIPPV interact with the spontaneous breathing of a VLBW infant (1).



Fig 1. The GIULIA Flow-sensor




Fig 2. Reading from the top, this recording shows the delivered pressure, flow and thoracic impedance of a VLBW infant treated with NIPPV. Note the interactions of NIPPV mandatory cycles (back-up rate of 20 breaths/min) with the spontaneous respiratory rhythm of the patient. The infant is not entrained with the ventilator and the mechanical cycles start (red lines) at different stages of the spontaneous breathing cycle: 1 peak of breath, 2 mid-expiration, 3 late expiration, 4 peak of breath, 5 early expiration, 6 mid-expiration.

Asynchronous mechanical breaths may induce laryngeal closure, alter spontaneous respiratory rhythm, increase WOB, increase abdominal distention, cause volutrauma and pneumothorax and have harmful effects on BP and CBF (2, 3).

Synchronous mechanical breaths may reduce inspiratory effort, increase ventilation, reduce breathing frequency, reduce thoraco-abdominal asynchrony and decrease WOB (3-7).

A common objection to the use of a flow sensor for non-invasive ventilation is that its reliability can be altered by the continuous flow passing through it due to the variable leaks from the infant's nostrils and mouth. To show that this objection is not valid, we used a simulated neonatal model to




Fig 3. This recording shows the interactions of flow-SNIPPV with spontaneous breathing of the newborn. The infant is now well entrained with the ventilator and mechanical breaths start (red lines) immediately after the beginning of the patient's spontaneous ones.


demonstrate both the reliability of our flow-sensor with different measured leaks through it and the performance of the GIULIA ventilator (8). The GIULIA flow-sensor detected 100% of the simulated spontaneous breaths in presence of any tested amount of leak from the prongs. The mean response time, measured from the beginning of inspiration to the beginning of the inspiratory pressure rise in the circuit, was  $64 \pm (SD) 7$  ms (range 46-77 ms). These data prove that the GIULIA flow-sensor can detect very small inspiratory volumes and flows, and that its performance is not affected by the amount of leaks. Another potential disadvantage of this device is the increase of dead volume, but this is only a theoretical problem since expiratory flow vents mainly from the patient's mouth.

# Giulia

Evidence from clinical trials indicates SNIPPV is more effective than nasal CPAP:

- in improving the success rate of extubation (9, 10)
- in reducing the need for intubation in the acute phase of RDS after surfactant (10-13)
- in infants with apnoea (1)

The GIULIA ventilator, initially developed only to carry out flow-SNIPPV, has recently been improved to perform also all the invasive ventilation modes (Fig 4).



## IMPROVEMENTS

The New GIULIA has been improved with:

- Invasive Ventilation Modes
- Encoder and color 10.4" Touch-Screen Digital Settings
- Manual and Automatic Alarms
- Pressure, Flow, Volume and Loop Graphical Trends
- Measurements of Compliance and Resistance
- Acoustic Signal of Trigger Activation
- O<sub>2</sub> Flush
- Battery Level

## CHARACTERISTICS

The New GIULIA can be used with patients of up to 5 kg of weight.

## VENTILATION MODES

The following conventional invasive and non-invasive ventilation modes can be set:

- CPAP - NCPAP
- SIMV - SNIMV
- SIPPV - SNIPPV
- IPPV - NIPPV

## FLOW SENSORS

GIULIA is supplied with two flow-sensors. Both are simple differential pressure transducers without any electrical components: one for invasive ventilation to calculate tidal volume, and the other for non-invasive ventilation, which comes in 2.5 and 3.5 mm sizes. The flow sensor for invasive ventilation guarantees a response time of less than 100 ms and is autoclavable and reusable.

## SMART FLOW KIT NIV

The new flow-sensor for non-invasive ventilation gives maximum lightness and comfort for the patient. The single-use "Smart-Flow NIV Kit" for non-invasive ventilation includes nasal prongs and bonnets, and comes in four different sizes for infants of different weights.

## NASAL PRONGS

The nasal prongs are made from a very soft material and are designed to adapt to the newborn's nostrils in the most comfortable way.

## ALARM SYSTEM

The GIULIA ventilator has all the alarms necessary for correct and safe management of invasive and non-invasive respiratory care. The alarms are both visual and acoustic, and are color-coded for priority.

## HUMIDIFIER

The GIULIA ventilator is compatible with all commercial humidifiers, however GINEVRI strongly recommends using it with the new WETTY humidifier, which ensures a high level of humidity in the respiratory circuit with a very low quantity of condensation.

## O<sub>2</sub> FLUSH

This function delivers a preset oxygen concentration for a predetermined time.

# Accessories

## WETTY

WETTY is a humidifier for heating and humidifying the gas in the patient's respiratory circuit. The gas temperature and the humidity are regulated by a servo-controlled system. The humidity can be set at five different levels.  
Code number 8049



## TROLLEY

This dedicated trolley allows GIULIA to be used easily with its accessories.  
Code number 11472A70

### OPTIONAL TROLLEY EXTRAS:

**DRAWERS.**  
Code number 11401A70  
**EXTRACTABLE SHELF.**  
Code number 11405A70  
**IV POLE.**  
Code number 6922



**TRACHEAL TEST-LUNG**  
Code number 11574A70



**NASAL PRONGS TEST-LUNG**  
Code number 11953A70



**TRACHEAL FLOW-SENSOR**  
(Autoclavable)  
Code number S103561300

# Consumables



## SMART FLOW KIT NIV (Disposable)

Four different sizes, color coded:

**RED KIT EXTRA-SMALL:** Nasal prongs Ø 2 mm, length 8 mm  
+ Flow-sensor, Ø 2.5 mm + Bonnet 25 cm. Code number 12898A08  
**GREEN KIT SMALL:** Nasal prongs Ø 2 mm, length 10 mm  
+ Flow-sensor, Ø 2.5 mm + Bonnet 25 cm. Code number 12898B08  
**WHITE KIT MEDIUM:** Nasal prongs Ø 3 mm, length 12 mm  
+ Flow-sensor, Ø 2.5 mm + Bonnet 30 cm. Code number 12898C08  
**BLUE KIT LARGE:** Nasal prongs Ø 4 mm, length 14 mm  
+ Flow-sensor, Ø 3.5 mm + Bonnet 35 cm. Code number 12898D08



**CONNECTION CIRCUIT FOR TRACHEAL FLOW-SENSOR**  
(Disposable)  
Code number 12936A08



**BONNETS**  
(Disposable)  
**RED EXTRA-SMALL:** 25 cm.  
Code number 11659A08  
**GREEN SMALL:** 25 cm.  
Code number 11659B08  
**WHITE MEDIUM:** 30 cm.  
Code number 11659C08  
**BLUE LARGE:** 35 cm.  
Code number 11659D08



**NASAL PRONGS**  
(Disposable)  
**EXTRA-SMALL:** Ø 2 mm, length 8 mm.  
Code number 12251A08  
**SMALL:** Ø 2 mm, length 10 mm.  
Code number 6968A08  
**MEDIUM:** Ø 3 mm, length 12 mm.  
Code number 6969A08  
**LARGE:** Ø 4 mm, length 14 mm.  
Code number 12205A08



**PATIENT CIRCUIT**  
(Disposable) Heated inspiratory pipes, with humidification chamber and water trap.  
Code number 12043A08



**O<sub>2</sub> SENSOR**  
Code number 10267A73



**EXPIRATORY VALVE MEMBRANE**  
Code number 11654A08



**DUST FILTER**  
Code number 12340A73

## References

- 1) Gizzi C, Montecchia F, Panetta V, Castellano C, Mariani C, Campelli M, Papoff P, Moretti C, Agostino R. "Is synchronized NIPPV more effective than NIPPV and NCPAP in treating Apnoea of Prematurity (AOP)? A randomised crossover trial", *Arch Dis Child Fetal Neonatal Ed* Published Online First: 15 Oct 2014 doi:10.1136/archdischild-2013-305892.
- 2) Greenough Z, Dimitriou G, Prendergast M, and Milner AD. "Synchronized mechanical ventilation for respiratory support in newborn infants", *Cochrane Database Syst Rev*, no. 1, Article ID CD000456, 2008
- 3) Chang HY, Claire N, D'Ugard C, Torres J, Nwajei P, and Bancalari E. "Effects of synchronization during nasal ventilation in clinically stable preterm infants", *Pediatric Research*, vol. 69, no. 1, pp. 84–89, 2011.
- 4) Kiciman NM, Andréasson B, Bernstein G, Mannino FL, Rich W, Henderson C, Heldt GP. "Thoracoabdominal motion in newborns during ventilation delivered by endotracheal tube or nasal prongs", *Pediatric Pulmonology*, vol. 25, no. 3, pp. 175–181, 1998.
- 5) Owen LS, Morley CJ, Dawson JA, and Davis PG. "Effects of non-synchronised nasal intermittent positive pressure ventilation on spontaneous breathing in preterm infants", *Archives of Disease in Childhood: Fetal and Neonatal Edition*, vol. 96, no. 6, pp. F422–F428, 2011.
- 6) Moretti C, Gizzi C, Papoff P, Lampariello S, Capoferri M, Calcagnini G, Bucci G. "Comparing the effects of nasal synchronized intermittent positive pressure ventilation (nSIPPV) and nasal continuous positive airway pressure (nCPAP) after extubation in very low birth weight infants", *Early Human Development*, vol. 56, no. 2-3, pp. 167–177, 1999.
- 7) Aghai ZH, Saslow JG, Nakhla T, Milcarek B, Hart J, Lawrysh-Plunkett R, Stahl G, Habib RH, Pyon KH. "Synchronized nasal intermittent positive pressure ventilation (SNIPPV) decreases work of breathing (WOB) in premature infants with respiratory distress syndrome (RDS) compared to nasal continuous positive airway pressure (NCPAP)", *Pediatric Pulmonology*, vol. 41, no. 9, pp. 875–881, 2006.
- 8) Moretti C, Papoff P, Gizzi C, Montecchia F, Giannini L, Fassi C, Midulla F, Agostino R, Sanchez-Luna M. "Flow-synchronized nasal intermittent positive pressure ventilation in the preterm infant: development of a project", *J Pediatr Neonat Individual Med*. 2013;2(2).
- 9) Moretti C, Giannini L, Fassi C, Gizzi C, Papoff P, Colarizi P. "Nasal flow-synchronized intermittent positive pressure ventilation to facilitate weaning in very low-birthweight infants: unmasked randomized controlled trial". *Pediatr Int* 2008; 50: 85–91.
- 10) Lemyre B, Davis PG, De Paoli AG, Kirpalani H. "Nasal intermittent positive pressure ventilation (NIPPV) versus nasal continuous positive airway pressure (NCPAP) for preterm neonates after extubation", *Cochrane Database Syst Rev*. 2014 Sep 4;9: CD003212. doi: 10.1002/14651858.CD003212.pub2.
- 11) Santin R, Brodsky N, Bhandari V. "A prospective observational pilot study of synchronized nasal intermittent positive pressure ventilation (SNIPPV) as a primary mode of ventilation in infants > or = 28 weeks with respiratory distress syndrome (RDS)", *J. Perinatol*. 2004; 24:487-93.
- 12) Bhandari V, Gavino RG, Nedrelow JH, Pallela P, Salvador A, Ehrenkranz RA, Brodsky NL. "A randomized controlled trial of synchronized nasal intermittent positive pressure ventilation in RDS", *J. Perinatol*, vol. 27, no. 11, pp. 697–703, 2007.
- 13) Gizzi C, Papoff P, Giordano I, Massenzi L, Barbàra CS, Campelli M, Panetta V, Agostino R, Moretti C. "Flow-synchronized nasal intermittent positive pressure ventilation for infants <32 weeks' gestation with respiratory distress syndrome". *Crit Care Res Pract*. 2012; 2012:301818.

## Technical Specifications

<b>General</b>	MPG–Class Dimensions (WxHxD) Weight	IIb 38x33x38 cm 18 kg		
<b>Power supply</b>	Mains Battery	230 V AC (50-60Hz),50VA 115 V AC (50- 60Hz),50VA Voltage	24VDC, 45min	
<b>Gas supply</b>		AIR O <sub>2</sub>	3–5 bar 3–5 bar	
<b>Operating mode</b>	Nasal Tracheal	With Nasal prongs With endotracheal tube		
<b>Ventilation modes</b> (INVASIVE & NON INVASIVE)		<b>SIMV, IMV ,SIPPV ,IPPV, PS, NCPAP,SNIPPV, NSIMV, NIPPV,BI-LEVEL</b>		
<b>Parameters</b>	Pmax PIP PEEP CPAP PSV Flow Inspiration time Expiration time FiO <sub>2</sub> Trigger T <sub>apn</sub> T <sub>bck</sub> Pre-oxygenation	0–100 cm H <sub>2</sub> O 10–100 cm H <sub>2</sub> O 0–50 cm H <sub>2</sub> O 0–15 cm H <sub>2</sub> O 0 – 80 cm H <sub>2</sub> O 1-40 L/min 0,1–5 s 0,4–60 s 21–100% Flow & Pressure 1–60s 1–180s 21-100 % FiO <sub>2</sub>		
<b>Monitoring</b>	<b>Display</b> <b>Pressure</b>  <b>Volume</b> V <sub>te</sub>  <b>Respiration frequency</b> <b>FiO<sub>2</sub></b> <b>Resistance</b> <b>Compliance</b> <b>Graph Display</b> <b>Data freeze and Loops</b>	10,4" Color–TFT PIP PEEP  0-500ml MinuteVolume  0–999 b/min 0–100% 0–9999cmH <sub>2</sub> O/l/s 0–999ml/cmH <sub>2</sub> O pressure , Volume ,Flow <b>P(t), V(t) , V'(t) , V(P) , V'(V), V'(P)</b>	-20 ± 100 cmH <sub>2</sub> O -20 ± 50 cmH <sub>2</sub> O   0-25 l/min   Accuracy: ± 5Vol	
<b>Alarms</b>	<b>Parameter</b> PIP PEEP V min VTE FiO <sub>2</sub>	<b>Lower Limit</b> From (PEEP+3) to (set PIP) From 0-50 or (PIP-2) 0.0– 6.5 0 – 20 20–100	<b>Upper limit</b> From (set PIP) to 100 0-50 0.0 – 25 0–500 22–100	<b>Unit</b> cmH <sub>2</sub> O cmH <sub>2</sub> O l/min ml %
<b>Data output</b>	USB			
<b>Flow Sensor</b>	Tracheal Pneumotachograph (Autoclavable / Reusable) Nasal Pneumotachograph (Reusable)	500g–5Kg Diam.2.5 Diam.3.5	500g–2Kg 2–5 Kg	



## Quality System

ISO 9001:2008  
ISO 13485:2012



The specifications in this catalogue are indicative. The company GINEVRI srl reserves the right to make changes, without further notice, to the products described within this catalogue in order to improve reliability, function or design.

**GINEVRI**

Quality for life

GINEVRI srl  
Via Cancelliera, 25/b  
00041 Albano Laziale  
(Rome) - Italy  
Tel.: +39 06 93459 330  
Fax: +39 06 93 459 393  
e-mail: [export@ginevri.com](mailto:export@ginevri.com)  
[www.ginevri.com](http://www.ginevri.com)

Authorized Distributor