

Aquatherm Aquatherm



GINEVRI
Quality for life

Aquatherm

A soft and warm nest for the little patient to make him feel like in the maternal womb...



When the baby is born, also the premature one, it requires different treatments, destinate to facilitate the growth phase and especially to avoid trauma that could affect the physical state of the baby.

At this stage it is very helpful to use the **Aquatherm** because it provides an additional beneficial source of heat to the little patient.

In addition, the **Baby Bundle**, flexible and easily shaped, is a gentle and effective aid to the newborn to find his natural position.

Technical specifications

Voltage	220 V \pm 10%
Frequency	50 Hz
Maximum Power	110 W
Current	0.5 A
Class / Type	Ia / B
Total weight	4.2 Kg
Dimension	
Thermo Driver	19 x 15 x 19 cm
Mattress	35 x 61 x 3 cm

Aquatherm

The standard configuration of the device **Aquatherm** has two units:



Thermo Driver

Unit that servo-controls the rate of heat supplied to the mattress in order it permits to achieve and maintain the temperature value to the little patient set by the operator. The **Thermo Driver** is provided with controls and alarms to ensure at the same time easy use and maximum reliability of the system **Aquatherm**.

Heating Mattress

The mattress unit is the result of a new technology "Aquagel" that offers all the advantages of a soft surface and impressed around the newborn, effective in minimizing the vibration and heat loss towards the surface and able to store the heat even when the device **Thermo Driver** is not powered or in case of internal transport. It is permeable to X-rays.



The **Aquatherm** system is portable, lightweight and easy positioning:

- **Thermo Driver** unit is equipped with legs that make it stable on any flat surface as well as a support to place it on vertical surfaces
- the mattress (heating unit and Aquagel) can be placed on a bench care (baby-changing table) as well as inside the incubator trays and however where you need a ground plane heated.

The **Aquatherm** system, by virtue of its property to store heat, it can be used in **Baby Crib Nido Ginevri** in the transport from a ward to another.



Since 1954 Ginevri is a leading designer, manufacturer and world-wide distributor of electromedical equipment for neonatal and pediatric care. Ginevri's policy has always been to promote the best quality, safety, user friendliness and easy maintenance of its products for the full satisfaction of its customers: public and private hospitals.

Quality System

ISO 9001
ISO 13485

The specifications in this catalogue are indicative. Ginevri, the company, reserves the right to make changes, without further notice, to the products described within this catalogue in order to improve reliability, function or design.

GINEVRI

Quality for life

GINEVRI srl
Via Cancelliera, 25/b
00041 - Cecchina di Albano
(Roma) - Italy

Tel.: +39 06 93 459 330
Fax: +39 06 93 459 393
e-mail: mail@ginevri.com
www.ginevri.com

Authorized Distributor

